

# Varilux Comfort® Family of Lenses



Shirley B. focusing  
on her balance.

**Varilux Comfort®** lenses are the gold standard for delivering superior vision by providing:

- Natural vision, Natural Posture
- Easy adaptation
- Ease of fitting

Shirley B. focusing  
on her students.

**Varilux® lenses**

**SUPERIOR VISION** WHATEVER THE DISTANCE.  
WHATEVER THE SITUATION.



**VARILUX**  
Comfort®

# Varilux Comfort® Family of Lenses



## VARILUX Comfort®

- Easy access to the near vision zone
- Wide fields of vision during daily activities
- Limited swim
- Minimum fitting height of 17mm



## VARILUX Comfort Short™

### All of the benefits of Varilux Comfort™ lenses plus:

- Wider distance vision than any competitive short corridor PAL
- Designed specifically for small frame wearers
- Fitting height of 14mm to 17mm



## VARILUX Comfort DRx™

## VARILUX Comfort Short DRx™

### All of the benefits of Varilux Comfort™ and Varilux Comfort Short™ lenses plus:

- Precisely places the near zone in the optimal position for patient's individual prescription
- Improved patient optics
- Essilor Digital Surfacing Process Control to ensure quality and performance
- Varilux Comfort DRx™ lenses - Minimum fitting height of 17mm
- Varilux Comfort Short DRx™ lenses - Fitting height of 14mm to 17mm



## VARILUX Comfort Enhanced™

### All of the benefits of Varilux Comfort™ and Varilux Comfort DRx™ lenses plus:

- Reduced peripheral distortion to expand visual fields throughout the lens
- Design elements on both sides of the lens through DualOptix™ technology
- Customized to your patient and the frame they choose
- Available in fitting heights from 14mm to 25mm

**Crizal®**

For the clearest vision possible, Varilux® lenses are recommended with Crizal® No-Glare lenses.

To learn more visit [VariluxUSA.com](http://VariluxUSA.com)

**VARILUX®**  
Natural Vision. Forever.